

SPENCER-WELLS HAEMOSTATIC FORCEPS

straight; serrated

20.0cm (8")



SKU: 040-290-20

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH338
KLS Martin	13-248-20
AUTHORS	SPENCER, WELLS
Medicon	15.76.20
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH338
KLS Martin	13-248-20

AUTHORS

[SPENCER](#), [WELLS](#)

Medicon

15.76.20

SPECIFICATIONS

[serrated](#), [straight](#)