

## SPENCER-WELLS HAEMOSTATIC FORCEPS

curved; serrated

14.0cm (5.6")




**SKU:** 040-291-14

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>KLS Martin</b>	13-249-14
<b>AUTHORS</b>	<a href="#">SPENCER</a> , <a href="#">WELLS</a>
<b>Medicon</b>	15.77.14
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>KLS Martin</b>	13-249-14
<b>AUTHORS</b>	<a href="#">SPENCER</a> , <a href="#">WELLS</a>
<b>Medicon</b>	15.77.14

## SPECIFICATIONS

[curved](#), [serrated](#)