

SPENCER-WELLS HAEMOSTATIC FORCEPS

curved; serrated

15.5cm (6.2")




SKU: 040-291-15

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH335
KLS Martin	13-249-15
AUTHORS	SPENCER, WELLS
Medicon	15.77.16
SPECIFICATIONS	curved , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH335
KLS Martin	13-249-15

AUTHORS

[SPENCER](#), [WELLS](#)

Medicon

15.77.16

SPECIFICATIONS

[curved](#), [serrated](#)