

## SPENCER-WELLS HAEMOSTATIC FORCEPS

curved; serrated

18.0cm (7.2")




**SKU:** 040-291-18

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	BH337
<b>KLS Martin</b>	13-249-18
<b>AUTHORS</b>	<a href="#">SPENCER, WELLS</a>
<b>Medicon</b>	15.77.18
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	BH337
<b>KLS Martin</b>	13-249-18

**AUTHORS**

[SPENCER](#), [WELLS](#)

**Medicon**

15.77.18

**SPECIFICATIONS**

[curved](#), [serrated](#)