

SPENCER-WELLS HAEMOSTATIC FORCEPS

straight; serrated

13.0cm (5.2")

SKU: 040-290-13

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH332
KLS Martin	13-248-13
AUTHORS	SPENCER, WELLS
Medicon	15.76.13
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH332
KLS Martin	13-248-13

AUTHORS

[SPENCER](#), [WELLS](#)

Medicon

15.76.13

SPECIFICATIONS

[serrated](#), [straight](#)