

## SPENCER-WELLS HAEMOSTATIC FORCEPS

straight; serrated

13.0cm (5.2")



**SKU:** 040-290-13

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                |
| <b>TYPE</b>           | <a href="#">Artery Forceps</a>                      |
| <b>Aesculap</b>       | BH332   |
| <b>KLS Martin</b>     | 13-248-13   |
| <b>AUTHORS</b>        | <a href="#">SPENCER</a> , <a href="#">WELLS</a>     |
| <b>Medicon</b>        | 15.76.13  |
| <b>SPECIFICATIONS</b> | <a href="#">serrated</a> , <a href="#">straight</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                   |                                |
|-------------------|--------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>           |
| <b>TYPE</b>       | <a href="#">Artery Forceps</a> |
| <b>Aesculap</b>   | BH332                          |
| <b>KLS Martin</b> | 13-248-13                      |

**AUTHORS**

[SPENCER](#), [WELLS](#)

**Medicon**

15.76.13

**SPECIFICATIONS**

[serrated](#), [straight](#)