

SPHERICAL BURR

5mm Borchardt shaft

Ø 4.0mm



SKU: 190-286-40

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-----------------------------|
| UNIT | pcs. |
| TYPE | Hand Drills |
| KLS Martin | 22-196-40 |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-----------------------------|
| UNIT | pcs. |
| TYPE | Hand Drills |
| KLS Martin | 22-196-40 |