

## SPHERICAL BURR; CONICAL

fine spacing; 5mm Borchardt shaft

Ø 3.0mm



**SKU:** 190-287-30

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hand Drills</a>
<b>KLS Martin</b>	22-198-30

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hand Drills</a>
<b>KLS Martin</b>	22-198-30