

SPLINTER FORCEPS

straight; serrated; without Pin

8.0cm (3.2")

SKU: 020-662-08

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Splinter Forceps
Aesculap	BD952
KLS Martin	12-604-08
Medicon	07.80.08
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Splinter Forceps
Aesculap	BD952
KLS Martin	12-604-08
Medicon	07.80.08

SPECIFICATIONS

[serrated](#), [straight](#)