

## SPLINTER FORCEPS

straight; serrated; without Pin

9.0cm (3.6")

**SKU:** 020-662-09

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Splinter Forceps</a>
<b>Aesculap</b>	BD301
<b>KLS Martin</b>	12-604-09
<b>Medicon</b>	07.80.09
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Splinter Forceps</a>
<b>Aesculap</b>	BD301
<b>KLS Martin</b>	12-604-09
<b>Medicon</b>	07.80.09

## SPECIFICATIONS

[serrated](#), [straight](#)