

SPLINTER FORCEPS

straight; serrated; with Pin

8.0cm (3,2")

SKU: 020-666-08

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Splinter Forceps](#)

SPECIFICATIONS [serrated](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Splinter Forceps](#)

SPECIFICATIONS [serrated](#), [straight](#)