

SPLINTER FORCEPS

straight; serrated; with Pin

10.5cm (4.2")

SKU: 020-666-10

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Splinter Forceps
Aesculap	BD302
KLS Martin	12-604-10
Medicon	07.80.10
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Splinter Forceps
Aesculap	BD302
KLS Martin	12-604-10
Medicon	07.80.10

SPECIFICATIONS

[serrated](#), [straight](#)