

SPLINTER FORCEPS

straight; serrated; with Pin

12.5cm (5")

SKU: 020-666-13

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Splinter Forceps](#)

Aesculap BD304

SPECIFICATIONS [serrated](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Splinter Forceps](#)

Aesculap BD304

SPECIFICATIONS [serrated](#), [straight](#)