

## SPLINTER FORCEPS

straight; serrated; with Pin

12.5cm (5")

**SKU:** 020-666-13

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Splinter Forceps](#)

**Aesculap** BD304

**SPECIFICATIONS** [serrated](#), [straight](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Splinter Forceps](#)

**Aesculap** BD304

**SPECIFICATIONS** [serrated](#), [straight](#)