

SPLINTER FORCEPS

curved; serrated; with Pin

8.0cm (3.2")

SKU: 020-667-08

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Splinter Forceps](#)

SPECIFICATIONS [curved](#), [serrated](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Splinter Forceps](#)

SPECIFICATIONS [curved](#), [serrated](#)