

## SPLINTER FORCEPS

curved; serrated; with Pin

10.5cm (4.2")

**SKU:** 020-667-10

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Splinter Forceps</a>
<b>Aesculap</b>	BD312
<b>KLS Martin</b>	12-605-10
<b>Medicon</b>	07.81.10
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Splinter Forceps</a>
<b>Aesculap</b>	BD312
<b>KLS Martin</b>	12-605-10
<b>Medicon</b>	07.81.10

## SPECIFICATIONS

[curved](#), [serrated](#)