

SPLINTER FORCEPS

curved; serrated; with Pin

13.0cm (5.2")

SKU: 020-667-13

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Splinter Forceps](#)

Aesculap BD314

SPECIFICATIONS [curved](#), [serrated](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Splinter Forceps](#)

Aesculap BD314

SPECIFICATIONS [curved](#), [serrated](#)