

SPLINTER FORCEPS

straight; smooth; without Pin

11.5cm (4.6")

SKU: 020-660-11

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Splinter Forceps](#)

KLS Martin 12-600-11

SPECIFICATIONS [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Splinter Forceps](#)

KLS Martin 12-600-11

SPECIFICATIONS [straight](#)