

## SPRAY-PROBE FOR COLOURING TECHNIQUE

round metal tip, type Spray (single use)

tubus 1.8 x 230cm, with mandrin



**SKU:** PK-8MREW2302230Q-S

**Category:** [Endoscopy](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)