

SPRAY-PROBE FOR COLOURING TECHNIQUE

round metal tip type Spray, double LL

1.8mm, length: 230cm

SKU: PK-8MRPK1802230MQ

Category: [Endoscopy](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Spray-Probes, Washing Tubes](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Spray-Probes, Washing Tubes](#)