

## SPRAY-PROBE FOR COLOURING TECHNIQUE

round metal tip type Spray, double LL

1.8mm, length: 230cm

**SKU:** PK-8MRPK1802230MQ

**Category:** [Endoscopy](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Spray-Probes, Washing Tubes](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Spray-Probes, Washing Tubes](#)