

STANDARD DRESSING FORCEPS

straight; serrated; narrow

16.0cm (6.4")

SKU: 020-004-16

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Dressing Forceps , Tissue Forceps |
| Aesculap | BD029 |
| KLS Martin | 12-120-16 |
| AUTHORS | STANDARD |
| Medicon | 06.01.16 |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|---|
| UNIT | pcs. |
| TYPE | Dressing Forceps , Tissue Forceps |
| Aesculap | BD029 |
| KLS Martin | 12-120-16 |

AUTHORS

[STANDARD](#)

Medicon

06.01.16

SPECIFICATIONS

[serrated](#), [straight](#)