

STANDARD DRESSING FORCEPS

straight; serrated; narrow

20.0cm (8")

SKU: 020-004-20

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD031
KLS Martin	12-120-20
AUTHORS	STANDARD
Medicon	06.01.20
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD031
KLS Martin	12-120-20

AUTHORS

[STANDARD](#)

Medicon

06.01.20

SPECIFICATIONS

[serrated](#), [straight](#)