

STANDARD DRESSING FORCEPS

straight; serrated; fine

11.5cm (4.6"); width of jaw 1.3 +/-0.2

SKU: 020-010-11

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
AUTHORS	STANDARD
Medicon	06.02.11
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
AUTHORS	STANDARD
Medicon	06.02.11
SPECIFICATIONS	serrated , straight