

## STANDARD DRESSING FORCEPS

straight; serrated; fine

11.5cm (4.6"); width of jaw 1.3 +/-0.2

**SKU:** 020-010-11

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>AUTHORS</b>	<a href="#">STANDARD</a>
<b>Medicon</b>	06.02.11
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>AUTHORS</b>	<a href="#">STANDARD</a>
<b>Medicon</b>	06.02.11
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>