

## STANDARD DRESSING FORCEPS

straight; serrated; fine

13.0cm (5.2"); width of jaw 1.3 +/-0.2

**SKU:** 020-010-13

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>KLS Martin</b>	12-126-13
<b>AUTHORS</b>	<a href="#">STANDARD</a>
<b>Medicon</b>	06.02.13
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>KLS Martin</b>	12-126-13
<b>AUTHORS</b>	<a href="#">STANDARD</a>
<b>Medicon</b>	06.02.13

## SPECIFICATIONS

[serrated](#), [straight](#)