

STANDARD DRESSING FORCEPS

straight; serrated; fine

18.0cm (7.2"); width of jaw 1.6 +/-0.2

SKU: 020-010-18

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD030
KLS Martin	12-120-18
AUTHORS	STANDARD
Medicon	06.02.18
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD030
KLS Martin	12-120-18

AUTHORS

[STANDARD](#)

Medicon

06.02.18

SPECIFICATIONS

[serrated](#), [straight](#)