

## STANDARD DRESSING FORCEPS

straight; serrated; fine

30.0cm (12"); width of jaw 2.5 +/-0.2

**SKU:** 020-010-30

**Categories:** [Forceps](#), [Open Surgery](#)



## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>AUTHORS</b>	<a href="#">MODELL USA</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>AUTHORS</b>	<a href="#">MODELL USA</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>