

## STANDARD DRESSING FORCEPS

straight; serrated

13.0cm (5.2"); width of jaw 2.7 +/- 0.2

**SKU:** 020-002-13

**Categories:** [Forceps](#), [Open Surgery](#)



## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>Aesculap</b>	BD045
<b>KLS Martin</b>	12-100-13
<b>AUTHORS</b>	<a href="#">STANDARD</a>
<b>Medicon</b>	06.00.13
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>Aesculap</b>	BD045
<b>KLS Martin</b>	12-100-13

## AUTHORS

[STANDARD](#)

## Medicon

06.00.13

## SPECIFICATIONS

[serrated](#), [straight](#)