

STANDARD DRESSING FORCEPS

straight; serrated

14.5cm (5.8"); width of jaw 3.1 +/-0.2

SKU: 020-002-14

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD047
KLS Martin	12-100-14
AUTHORS	STANDARD
Medicon	06.00.14
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD047
KLS Martin	12-100-14

AUTHORS

[STANDARD](#)

Medicon

06.00.14

SPECIFICATIONS

[serrated](#), [straight](#)