

STANDARD DRESSING FORCEPS

straight; serrated

16.0cm (6.4"); width of jaw 3.7 +/-0.2

SKU: 020-002-16

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD049
KLS Martin	12-100-16
AUTHORS	STANDARD
Medicon	06.00.16
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD049
KLS Martin	12-100-16

AUTHORS

[STANDARD](#)

Medicon

06.00.16

SPECIFICATIONS

[serrated](#), [straight](#)