

STANDARD DRESSING FORCEPS

straight; serrated

25.0cm (10"); width of jaw 3.7 +/-0.2

SKU: 020-002-25

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD052
KLS Martin	12-100-25
AUTHORS	STANDARD
Medicon	06.00.25
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD052
KLS Martin	12-100-25

AUTHORS

[STANDARD](#)

Medicon

06.00.25

SPECIFICATIONS

[serrated](#), [straight](#)