

## STANDARD DRESSING FORCEPS

straight; serrated

30.0cm (12"); width of jaw 5.0 +/-0.2

**SKU:** 020-002-30

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>Aesculap</b>	BD053
<b>KLS Martin</b>	12-100-30
<b>AUTHORS</b>	<a href="#">STANDARD</a>
<b>Medicon</b>	06.00.30
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>Aesculap</b>	BD053
<b>KLS Martin</b>	12-100-30

## AUTHORS

[STANDARD](#)

## Medicon

06.00.30

## SPECIFICATIONS

[serrated](#), [straight](#)