

## STANDARD FIBER COLD LIGHT GUIDE CABLE

without connectors

Ø 3,5 x 3000 mm Length

Tube colour grey; careful with Xenon light

**STANDARD**  
COLD LIGHT LINE



**SKU:** 540-803-30

**Categories:** [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Light Guide Cables</a>
<b>TYPE</b>	<a href="#">Standard (Cold Light)</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Light Guide Cables</a>
<b>TYPE</b>	<a href="#">Standard (Cold Light)</a>