

## STANDARD FIBER COLD LIGHT GUIDE CABLE

without connectors

Ø 4,9 x 2300 mm Length, Tube colour grey

careful with Xenon lights

**STANDARD**  
COLD LIGHT LINE



**SKU:** 540-805-23

**Categories:** [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Light Guide Cables</a>
<b>TYPE</b>	<a href="#">Standard (Cold Light)</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Light Guide Cables</a>
<b>TYPE</b>	<a href="#">Standard (Cold Light)</a>