

STANDARD FIBER COLD LIGHT GUIDE CABLE

without connectors

Ø 3,5 x 1800 mm Length

Tube colour grey; careful with Xenon light

STANDARD
COLD LIGHT LINE



SKU: 540-803-18

Categories: [Arthroscopy](#), [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	Standard (Cold Light)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Light Guide Cables
TYPE	Standard (Cold Light)