

STANDARD SHEATH;17FR.; ROTABLE; 1 STOPCOC

for bipolar working element 582-180-00

QL connection, colour code green



SKU: 582-180-15

Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|-----------------------------------|
| UNIT | pcs. |
| RANGE | Mini Resectoscopy |
| TYPE | Sheaths |

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|-----------------------------------|
| UNIT | pcs. |
| RANGE | Mini Resectoscopy |
| TYPE | Sheaths |