

## STANDARD TISSUE FORCEPS

straight; 2X3 teeth

30.0cm (12"); width of jaw 4.0 +/-0.2

**SKU:** 020-252-30

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>KLS Martin</b>	12-302-30
<b>AUTHORS</b>	<a href="#">STANDARD</a>
<b>SPECIFICATIONS</b>	<a href="#">straight</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>KLS Martin</b>	12-302-30
<b>AUTHORS</b>	<a href="#">STANDARD</a>
<b>SPECIFICATIONS</b>	<a href="#">straight</a> , <a href="#">toothed</a>