

STEINMANN EXTENSION PIN

with trocar point; triangular shaft

Ø 4.0mm; 16.0cm (6.4")



SKU: 212-348-16

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Steinmann pins
Aesculap	KX116
KLS Martin	22-562-16
AUTHORS	STEINMANN

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Steinmann pins
Aesculap	KX116
KLS Martin	22-562-16
AUTHORS	STEINMANN