

STEINMANN EXTENSION PIN

with trocar point; triangular shaft

Ø 4.0mm; 18.0cm (7.2")



SKU: 212-348-18

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Steinmann pins
Aesculap	KX118
KLS Martin	22-562-18
AUTHORS	STEINMANN

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Steinmann pins
Aesculap	KX118
KLS Martin	22-562-18
AUTHORS	STEINMANN