

## STEINMANN EXTENSION PIN

with trocar point; triangular shaft

Ø 4.0mm; 14.0cm (5.6")



**SKU:** 212-348-14

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Steinmann pins</a>
<b>Aesculap</b>	KX114
<b>KLS Martin</b>	22-562-14
<b>AUTHORS</b>	<a href="#">STEINMANN</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Steinmann pins</a>
<b>Aesculap</b>	KX114
<b>KLS Martin</b>	22-562-14
<b>AUTHORS</b>	<a href="#">STEINMANN</a>