

STEINMANN TENDON HOLDING FORCEPS

straight; serrated

15.5cm (6.2")




SKU: 195-246-16

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Seizing Instruments
KLS Martin	23-674-16
AUTHORS	STEINMANN
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Seizing Instruments
KLS Martin	23-674-16
AUTHORS	STEINMANN
SPECIFICATIONS	serrated , straight