

## STEINMANN TENDON HOLDING FORCEPS

straight; serrated

15.5cm (6.2")




**SKU:** 195-246-16

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tendon Seizing Instruments</a>
<b>KLS Martin</b>	23-674-16
<b>AUTHORS</b>	<a href="#">STEINMANN</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tendon Seizing Instruments</a>
<b>KLS Martin</b>	23-674-16
<b>AUTHORS</b>	<a href="#">STEINMANN</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>