

## STILLE DRESSING FORCEPS

straight; serrated

15.0cm (6"); width of jaw 2.7 +/-0.2

**SKU:** 020-008-15

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>Aesculap</b>	BD067
<b>KLS Martin</b>	12-130-15
<b>AUTHORS</b>	<a href="#">STANDARD</a>
<b>Medicon</b>	06.44.15
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>Aesculap</b>	BD067
<b>KLS Martin</b>	12-130-15

## AUTHORS

[STANDARD](#)

## Medicon

06.44.15

## SPECIFICATIONS

[serrated](#), [straight](#)