

## STILLE-LUER BONE RONGEUR

straight; wide jaw; double convertible

27.0cm (10.8")


**SKU:** 196-354-27

**Categories:** [Open Surgery](#), [Orthopaedics](#)



### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                          |
| <b>TYPE</b>           | <a href="#">Bone Rongeurs</a>                 |
| <b>KLS Martin</b>     | 23-912-27                                     |
| <b>AUTHORS</b>        | <a href="#">LUER</a> , <a href="#">STILLE</a> |
| <b>SPECIFICATIONS</b> | <a href="#">straight</a>                      |



### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                          |
| <b>TYPE</b>           | <a href="#">Bone Rongeurs</a>                 |
| <b>KLS Martin</b>     | 23-912-27                                     |
| <b>AUTHORS</b>        | <a href="#">LUER</a> , <a href="#">STILLE</a> |
| <b>SPECIFICATIONS</b> | <a href="#">straight</a>                      |