


## STOCKMANN PENIS CLAMP

8.0cm (3.2")



**SKU:** 170-345-08

**Categories:** [Open Surgery](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Foreign Body Forceps, Penis Clamps</a>
<b>Aesculap</b>	EF400
<b>KLS Martin</b>	31-875-08
<b>AUTHORS</b>	<a href="#">STOCKMANN</a>
<b>Medicon</b>	51.75.05

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Foreign Body Forceps, Penis Clamps</a>
<b>Aesculap</b>	EF400
<b>KLS Martin</b>	31-875-08
<b>AUTHORS</b>	<a href="#">STOCKMANN</a>
<b>Medicon</b>	51.75.05