

## STORAGE AND STERILISATION RACK

for 10 laparoscopy instruments Ø 3.0mm

Ø 5.0mm or Ø 10.0mm; 525x245x175mm

incl. 2 silicon holders; for container 580x280x200mm



**SKU:** 499-101-20

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Accessories](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Accessories](#)