

## STORAGE AND STERILISATION RACK

for 15 laparoscopy instruments with

Ø5.0mm or Ø10.0mm; 470x250x150mm



**SKU:** 499-100-25

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Accessories](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Accessories](#)