

## STRUEMPEL EAR FORCEPS

fenestrated

shaft 8.5cm (3.4")



**SKU:** 111-292-08

**Categories:** [Open Surgery](#), [Otology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Ear Forceps</a>
<b>Aesculap</b>	OF430
<b>KLS Martin</b>	36-448-08
<b>AUTHORS</b>	<a href="#">STRUEMPEL</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Ear Forceps</a>
<b>Aesculap</b>	OF430
<b>KLS Martin</b>	36-448-08
<b>AUTHORS</b>	<a href="#">STRUEMPEL</a>