

SUBRAMANIAN PEDIATRIC VASCULAR CLAMP

right angled; graduated; atraumatic

16.0cm (6.2")




SKU: 142-631-01

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
KLS Martin	24-561-16
AUTHORS	SUBRAMANIAN
SPECIFICATIONS	angled



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
KLS Martin	24-561-16
AUTHORS	SUBRAMANIAN
SPECIFICATIONS	angled